ACRYL HOT (GB)

Material: ACRYL HOT may be processed by all polymerization techniques, but has been developed especially for the rapid polymerization by packing and injection procedure.

Mixing ratio: Recommended ratio: 2.5 parts powder (by weight): 1 part liquid (by weight).

Working time: About 1 hour, depending upon mixing ratio.

Working procedure: You may mix powder and liquid without following exact ratios. Add as much powder to liquid as is immediately absorbed. Stir thoroughly with wide spatula for 30 sec. After 10 min. the dough is workable for around 45 min. If longer workability is desired we recommend the following ratios be used:

- 2.0 part powder (by weight)

- 1.0 part liquid (by weight)

Doughed acrylic can be kept for around 24 hours in a refrigerator. At normal room temperature the working time will be around 60 min., at ratio 2.5 : 1 and 2 - 3 hours at ratio 2.0 : 1. The mechanical properties do not change.

Polymerization: ACRYL HOT can be cured by normal and long polymerization. Put flask into boiling water and polymerize for 20 minutes. Thicker sectioned dentures require that flask be placed in boiling water. Switch off heating. Wait for 15 minutes. Then polymerize for a further 10 minutes. Cool down flask as usual.

Important: The above time are based at room temperature of 23°C: higher temperatures will shorten the time lower temperature increases them.

Working with injection units: Process the flasks and plaster surfaces as described in the instructions for this system. ACRYL HOT is mixed as described for the injection-moulding technique, placed in the cylinder and injected. The polymerisation cycle is the same as that used for injection-moulding.

Repairs: For Repairs use ACRYL COLD. Before applying the acrylic dough, roughen and clean the fracture surfaces and coat them with monomer.

Storage: At room temperature, approx. 21° C. Please observe the shelf-life.

Shade: Pink-transparent.

Presentation: 500g and 1000g Powder or 500ml and 1000ml Liquid.

Side-effects: With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please contact us.

Contra-indications/interactions: We will supply the composition of our medical device upon request. If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross reactions of the product with other materials already in the patient's mouth before using the product.

Note: Please supply the dentist with the above information, if this medical device is used to produce a special model. Please consult the safety data sheets.

TROUBLESHOOTING		
PROBLEM	CAUSE	REMEDY
Acrylic does not fully cure	Incorrect mixing ratio	Adhere to the instructions for use Adhere to the mixing ratio
Acrylic is whitish on the fitting surface	Incorrect or poor application of separating agent	Apply ample separating agent Check separating agent if necessary Soak the model for approx. 10 min. in hand-hot water
Acrylic form bubbles on the fitting surface	Model not soaked	Soak the model for approx. 10 min. in hand-hot water
Acrylic surface white/porous	Model not soaked, recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Colour differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic/denture has whitish streaks	Incorrectly mixed/too much powder used	Thoroughly spatulate the acrylic dough/adhere to the mixing ratio

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Henry Schein U.K. Holdings Ltd.